

Towards precise walking strides segmentation from accelerometry data in free-living: ADEPT method and further developments

Marta Karas

Advised by Ciprian Crainiceanu, Jacek Urbanek

BIRS: Use of Wearable and Implantable Devices in Health Research
Feb 25, 2020

Wearable and Implantable Technology (WIT) group



[Ciprian Crainiceanu](#),
[Jacek Urbanek](#)




[Vadim Zipunnikov](#)



[Andrew Leroux](#)




[Lacey Etzkorn](#)



[Marta Karas](#)



[Debangnan Dey](#)



[Erjia Cui](#)



[Jennifer Xu](#)



[Jawei Bai](#)




[John Muschelli](#)





Outline

1. ADEPT method for pattern segmentation, with application to walking strides




2. Towards walking segmentation in free-living: digital biomarkers of pain/stiffness in arthritis patients




Raw accelerometry data

Macro-scale – summarized data (e.g. 1 min intervals)



- Activity counts, number of steps etc.



Micro-scale – raw accelerometry data (10Hz+)




Raw accelerometry data of walking



ADEPT method for automatic pattern segmentation from acc. time-series

Adaptive Empirical Pattern Transformation (**ADEPT**) method optimized to segment walking strides from acc vector magnitude time-series:

1. Use **stride pattern template**:
 - scale (stretch/shrink),
 - translate (move)
2. Identify scale and time where template has the **highest similarity** (correlation, covariance) to accelerometry time-series




$$W_{\Psi}(s, \tau) = \int_{-\infty}^{\infty} x(t) \frac{1}{\sqrt{s}} \Psi\left(\frac{t-\tau}{s}\right) dt.$$

Template: ADEPT



Accelerometry
vector magnitude



ADEPT method for automatic segmentation of walking strides: strengths

- Uses data [r] vector magnitude for pattern-matching

⇒ **Method is invariant to sensor rotation (does not rely on device axes orientation relative to human body)**

- In pattern matching, considers simultaneously: (a) various pattern width values, (b) multiple stride patterns,

⇒ **Accounts for that different people walk differently and with various pace**

- Uses maxima-fine tuning step in the algorithm

⇒ **Algorithm returns precise location of start/end of a walking stride**

ADEPT method: application examples

Context:



- Supervised walking, semi-supervised walking, free-living walking [*work in progress*]
- Clinical intervention, clinical trial, epidemiological study

Derive measurements of walking/gait:

- **Cadence/walking stride duration**
- Functional data analysis of subject-specific **stride pattern**
- Diurnal patterns of walking
- Variability of those parameters within a trial

Use to:

- Characterize/differentiate individuals in population
- Characterize *changes* within those parameters for individual over time



Adaptive empirical pattern transformation (ADEPT) with application to walking stride segmentation

Marta Karas ✉, Marcin Strączkiewicz, William Fadel, Jaroslaw Harezlak, Ciprian M Crainiceanu, Jacek K Urbanek

Biostatistics, kxz033, <https://doi.org/10.1093/biostatistics/kxz033>

Published: 23 September 2019 **Article history** ▼

- Method application to segment walking strides from continuous walking from **32 healthy adults at 4 locations**: wrist, hip and both ankles
- Highest grade of results reproducibility (R)
 - All data and code available, results 100% reproducible

FEBRUARY 2020 WEB VERSION | ACTIGRAPHCORP.COM

RESEARCH SPOTLIGHT

BY ACTIGRAPH



[Adaptive Empirical Pattern Transformation \(ADEPT\) with Application to Walking Stride Segmentation](#)

Quantifying gait parameters and ambulatory monitoring of changes in these parameters have become increasingly important in epidemiological and clinical studies. Using high-density accelerometry measurements...[Read More](#)

ADEPT method: software



R package *adept* implements ADEPT method

- [CRAN index](#), [GitHub repo](#)
- Vignette 1: [Introduction to adept package](#)
- Vignette 2: [Walking strides segmentation with adept](#)
- Selected in Top new R packages in May 2019

Also:

- Language-agnostic pseudocode provided [[link](#)]
- **Implemented on GPU** for superior speed performance

May 2019: "Top 40" New CRAN Packages


📅 2019-06-25

by Joseph Rickert

Two hundred twenty-two new packages made it to CRAN in May, and it was more of an effort than usual to select the "Top 40". Nevertheless, here they are in nine categories, Computational Methods, Data, Machine Learning, Mathematics, Medicine, Science, Statistics, Utilities and Visualization.


Medicine

[adept v1.1.2](#): Provides functions for analyzing high-density data from walking strides collected from a wearable accelerometer worn during continuous walking activity. There is an [Introduction](#) and a [vignette](#) on walking stride segmentation.



Digital biomarkers of pain/stiffness in arthritis patients from wrist-collected accelerometry data

- **Arthritis** - swelling and tenderness of one or more of your joints; main symptoms are **joint pain and stiffness**, traditionally measured with questionnaires
- Issue: in drug development, we want to have pain and stiffness being **objectively measured**



© MAYO FOUNDATION FOR MEDICAL EDUCATION AND RESEARCH. ALL RIGHTS RESERVED.

Image source:

<https://www.mayoclinic.org/diseases-conditions/arthritis/symptoms-causes/>

Digital biomarkers of pain/stiffness in arthritis patients from wrist-collected accelerometry data

- **Arthritis** - swelling and tenderness of one or more of your joints; main symptoms are **joint pain and stiffness**, traditionally measured with questionnaires
- Issue: in drug development, we want to have pain and stiffness being **objectively measured**

Previous work (observational study):

- **5x sit-to-stand exercise time is associated with self-reported pain and stiffness scores**
 - Used wrist-worn sensor raw data to measure 5x sit-to-stand time



Digital biomarkers of pain/stiffness in arthritis patients from wrist-collected accelerometry data

- **Arthritis** - swelling and tenderness of one or more of your joints; main symptoms are **joint pain and stiffness**, traditionally measured with questionnaires
- Issue: in drug development, we want to have pain and stiffness being **objectively measured**

Previous work (observational study):

- **5x sit-to-stand exercise time is associated with self-reported pain and stiffness scores**
 - Used wrist-worn sensor raw data to measure 5x sit-to-stand time

We propose:

- **Segment walking strides from accelerometry data collected in free-living**
- Use data features of walking strides (stride duration time, others) as pain/stiffness digital biomarker candidates






Segment walking in free-living: algorithm

1. Use ADEPT pattern-matching algorithm to **identify individual walking strides** pattern

2. Identify **valid walking bouts** in free-living data


3. Derive data features for each walking stride data segment



Identified walking stride
(2 subsequent steps)


Walking bout:

- sequence of ≥ 3 walking **strides**
- subsequent **strides** must be “valid and similar to each other”



Describe each walking stride:


1. Duration time [s]
2. Vector magnitude peak-to-peak
3. Vector magnitude standard deviation
4. Vector magnitude count
5. Azimuth peak-to-peak
6. Azimuth standard deviation



Digital biomarkers of pain/stiffness in arthritis patients from wrist-collected accelerometry data: challenges


- **Validation of method for walking segmentation** from data collected at wrist in free-living:
 - No gold standard
 - No free-living raw accelerometry data with precise walking labels
- **Statistical modeling** of association between walking data features and self-reported scores of pain and stiffness
 - Account for context of walking?
 - Use walking information from morning only / whole day?
 - Use 1-point summary of per-day per-subject measurements, or use functional predictor?

Idea #1: Example: one study participant.



Idea #2:

- Right upper figure: based on other Novartis study, data collected with accelerometry at **belt buckle-worn sensor**.
- Right lower figure: Presented study, data collected with accelerometry non-dominant **wrist-worn sensor**.



Thank you!

